**T1**

**1000 mg/L:** Na<sup>+</sup>, K<sup>+</sup>, SO<sub>4</sub><sup>2-</sup>, Ca<sup>2+</sup>, Mg<sup>2+</sup>

**500 mg/L:** NO<sub>3</sub><sup>-</sup>

**200 mg/L:** Cl<sup>-</sup>

**100 mg/L:** PO<sub>4</sub><sup>3-</sup>

**50 mg/L:** Fe<sup>2+</sup>, Mn<sup>2+</sup>

**10 mg/L:** Fe<sup>3+</sup>

**1 mg/L:** NO<sub>2</sub><sup>-</sup>, Cu<sup>2+</sup>

### Datatablel / Data table

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**GB**

### LCK 330 Molybdenum

**Please check the "Edition Date" (see data table).**

**Safety advice and expiry date on package.**

#### Principle

Molybdate ions react with thioglycolic acid to form a yellow coloured complex, which is then photometrically evaluated.

#### Range of Application

Raw and ground water, drinking water, boiler water, process analysis

#### Interferences

The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions.

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

#### Removal of Interferences

The sample should be colourless and free of turbidities. Slight colourations can be taken into account with the help of a samplespecific blank reading. (see procedure **P1**)

If such colouration is due to organic impurities, these can be eliminated by digesting the sample with the Crack Set (LCW 902).

Turbidities are eliminated by filtration through a membrane filter (LCW 904).

**pH** sample .....4 – 9  
**Temperature** sample/reagents .....22°C  
**In case of not working at the right recommended temperature an incorrect result may be obtained.**

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
**Specific blank with coloured sample**

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
Pipette into a blank-value cuvette (LCW 919) **2.0 mL** distilled water.

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


Kuvet sluiten en zwenken.

Close cuvette and invert a few times.

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After **10 min** thoroughly clean the outside of the cuvette and evaluate.

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If the sample is coloured the blank-value must be measured and subtracted from the value obtained in the analysis (see procedure **P1**). Please observe the measuring range limits.

<sup>1)</sup> LASA 50 / 100  
XION 500  
CADAS 30 / 50 / 30S / 50S / 200 Barcode  
ISIS 9000  
DR 2800 / DR 3800 / DR 3900 / DR 5000 / DR 6000

	<b>Analysenküvette ①</b> <b>Cuve d'analyse</b> <b>Cuvetta d'analisi</b> <b>Analyse-kuvet</b> <b>Sample cuvette</b>
	✓

	Filter ① Filtre Filtro Filter Filter	Eprom ②	Test ③ - anwählen - choisir - selezionare - oproepen - select	Kontrollnr. ④ No. de contrôle No. di controllo Controlegetal Control no.	Analysenküvette ⑤ Cuve d'analyse Cuvetta d'analisi Analyse-kuvet Sample cuvette
<b>LASA 1 / plus</b>	330 nm	--	Mo LCK 330	1	✓
<b>LASA 10 / 20</b>	--	_ : 48	Mo LCK 330	1	✓

	Filter ① Filtre Filtro Filter Filter	Test ② - anwählen - choisir - selezionare - oproepen - select	Kontrollnr. ③ No. de contrôle No. di controllo Controlegetal Control no.	Leerwert (Probe) ④ Valeur à blanc (échantillon) Bianco (campione) Blanko (monster) Blank-value (sample)	Analysenküvette ⑤ Cuve d'analyse Cuvetta d'analisi Analyse-kuvet Sample cuvette
<b>LP2W</b>	340 nm	Test _ _	1	 	 
				LCW 919	✓

	Filter ① Filtre Filtro Filter Filter	Eprom ②	Mode ③ 	Test ④ - anwählen - choisir - selezionare - oproepen - select	Kontrollnr. ⑤ No. de contrôle No. di controllo Controlegetal Control no.	Analysenküvette, grüne Taste / Messen ⑥ Cuve d'analyse, touche verte / Mesurer Cuvetta d'analisi, tasto verde / Lettura Analyse-kuvet, groene toets / Meten Sample cuvette, green key / Read
<b>CADAS 200 Basis</b>	--	_ : 48	--	330	1	✓
<b>ISIS 6000</b>	--	_ : 48	<sup>2)</sup>	330	1	✓
<b>LASA 30</b>	340 nm	--	Dr. Lange	330	1	✓
<b>DR 1900</b>	--	--	<sup>3)</sup>	330	1	✓

<sup>2)</sup> KÜVETTEN-TEST

<sup>3)</sup> BARCODE-PROGRAMME

<sup>2)</sup> TEST EN CUVE

<sup>3)</sup> PROGR. CODE BARRE

<sup>2)</sup> CUVETTE-TEST

<sup>3)</sup> PROGRAMMI COD. A BARRE

<sup>2)</sup> KUVETTENTEST

<sup>3)</sup> BARCODEPROGRAMMA'S

<sup>2)</sup> CUVETTE TEST

<sup>3)</sup> BARCODE PROGRAMS

	Mode ①	Symbol ② Symbole Simbolo Symbool Symbol	Kontrollnr. ③ No. de contrôle No. di controllo Controlegetal Control no.	Leerwert (Probe) ④ Valeur à blanc (échantillon) Bianco (campione) Blanko (monster) Blank-value (sample)	Analysenküvette ⑤ Cuve d'analyse Cuvetta d'analisi Analyse-kuvet Sample cuvette
<b>CADAS 100 LPG210</b>	TEST	330	1	 	 
				LCW 919	✓