T1

1700 mg/L: Cl ⁻
1500 mg/L: Na ⁺ , K ⁺ , NH ₄ ⁺ , Ca ²⁺
300 mg/L: Cu ²⁺
200 mg/L: SO ₄ ²⁻
100 mg/L: NO ₃ ⁻
50 mg/L: SO ₃ ²⁻ , NO ₂ ⁻ , CN ⁻ , J ⁻ , CH ₃ COO ⁻ , Al ³⁺ , Pb ²⁺ , Mn ²⁺ , Cr ⁶⁺ , Sn ²⁺ , Fe ²⁺ , Zn ²⁺ , Hg ²⁺ , Cd ²⁺ , Ni ²⁺ , Ag ⁺
20 mg/L: Co ²⁺
10 mg/L: Fe ³⁺ , S ²⁻ , CH ₂ O
2 mg/L: Cr ³⁺
1 mg/L: H ₂ O ₂
5 Vol. %: CH ₃ OH, C ₂ H ₅ OH, (CH ₃) ₂ CO

Datatablel / Data table

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*) **Fenol I – Meetbereik I**
Phenols I – Measuring range I

) **Fenol II – Meetbereik II
Phenols II – Measuring range II

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GB

LCK 346 Phenols

! **Please check the "Edition Date" (see data table).**
Safety advice and expiry date on package.

Principle

In the presence of an oxidizing agent ortho- and meta-substituted phenols form coloured complexes with 4-aminoantipyrine (AAP).

Range of Application

Determination of the total phenols content of waste water, exhaust air (after absorption), and exhaust air condensates formed during the manufacture and processing of benzene, petroleum products, glass and mineral fibres, hardboard, coke, oil shale, hazardous waste, town gas, coal and brown coal products, tar, asphalt and bitumen.

Measuring range

	LASA aqua LASA 1/plus	LASA 20 LP1W/LP2W	LASA 30 Barcode, ISIS 6000 CADAS 200 Basis	CADAS 100
Phenols, Meas. range I	5 – 20 mg/L	5 – 50 mg/L	5 – 40 mg/L	5 – 30 mg/L
Phenols, Meas. range II	20 – 60 mg/L	20 – 200 mg/L	20 – 150 mg/L	20 – 100 mg/L

Interferences

The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions. Larger amounts of cobalt, iron(III), chromium(III), and sulphide cause high-bias results. Higher volume percentages of water-soluble organic solvents can cause high or low-bias results with differently substituted phenols. High concentrations of strong oxidizing and reducing agents in the sample interfere with the reaction process and must be removed before the analysis should be carried out. Other substances capable of combining with 4-aminoantipyrine (e.g. naphthols and aromatic amines) are also partially analysed and can simulate a higher phenol concentration. The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

Special note

If the measurements cannot be carried out immediately, the phenol samples should be cooled until they can be analysed.

pH sample2 – 11
Strongly acidic or alkaline samples must previously be adjusted to the correct pH.
Temperature sample/reagents20 – 25°C

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