**T1**

**1000 mg/l:** K<sup>+</sup>, Na<sup>+</sup>, Cl<sup>-</sup>

**500 mg/l:** NH<sub>4</sub><sup>+</sup>, SO<sub>4</sub><sup>2-</sup>, NO<sub>3</sub><sup>-</sup>, NO<sub>2</sub><sup>-</sup>, Mg<sup>2+</sup>

**200 mg/l:** Cu<sup>2+</sup>, Ni<sup>2+</sup>, Zn<sup>2+</sup>

**100 mg/l:** Ca<sup>2+</sup>, Al<sup>3+</sup>

**T2**

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**GB**

### LCK 433 Nonionic surfactants

**Please check the "Edition Date" (see data table).**

**Safety advice and expiry date on package.**

#### Principle

Nonionic surfactants (ethoxylates with 3 – 20 ether bridges) react with the indicator TBPE to form complexes, which are extracted in dichloromethane and photometrically evaluated.

#### Range of Application

Water, waste water, surface water, formulations, degreasing baths, wash solutions, process analysis

#### Interferences

Cationic surfactants cause high-bias results to be obtained.  
Anionic surfactants cause low-bias results to be obtained, depending on the concentration of the anionic surfactant. At anionic surfactant concentrations of 2.0 mg/l the bias is approx. 10%, while at concentrations of 20 mg/l the bias is between 40%.  
APGs (alkyl polyglycosides) are not measured.  
The ions listed in **T1** have been individually checked up to the given concentrations and do not cause interference. We have not determined cumulative effects and the influence of other ions.  
The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

#### For special attention

This method of surfactant analysis determines the sum of the surfactant types alkylphenol ethoxylates (AP(EO)<sub>n</sub>), fatty alcohol ethoxylates (FA(EO)<sub>n</sub>) and polyethylene glycols (PEG). Unless the composition of the sample is known, it is not possible to determine the individual substances. The result of the measurement is therefore expressed as a concentration of a standard surfactant with 10 ethoxy bridges as **mg/l TRITON x 100**.  
The result can only be converted to concentrations of other surfactants with the help of a special calibration (see **T2**).

**pH** sample ..... 4 – 9  
**Temperature** sample/reagents ..... 20 – 23°C  
**In case of not working at the right recommended temperature an incorrect result may be obtained.**

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
**Attention:**

The phase separation may take several hours if the cuvette is inverted or shaken too vigorously!

The phase separation must be complete before the Cuvette Test is evaluated. Incomplete phase separation causes high-bias results to be obtained.

If the cuvette is inverted too cautiously the extraction may be incomplete, so that low-bias results are obtained.

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
Close the cuvette and hold it between finger and thumb. **Mix intensively (2 - 3 times per second) for 2 min - do not shake.**

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
Then allow the cuvette to stand. Wait for the phase separation to finish (**approx. 2 min**).

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Thoroughly clean the outside of the cuvette and evaluate.

	<b>Analysenküvette ①</b> <b>Cuve d'analyse</b> <b>Cuvetta d'analisi</b> <b>Analyse-kuvet</b> <b>Sample cuvette</b>
 Barcode <sup>1)</sup>	✓

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