T1

NL

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EN

LCK 514 COD Chemical Oxygen Demand

! **Please check the "Edition Date" (see data table) and read the "Note". Safety advice and expiry date on package.**

Principe

Oxidizable substances react with sulphuric acid – potassium dichromate solution in the presence of silver sulphate as a catalyst. Chloride is masked by mercury sulphate. The green coloration of Cr³⁺ is evaluated.

Range of Application

Waste water, process analysis

Interferences

The method can be used for samples (or diluted samples) with chloride concentrations of up to 1500 mg/L. The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

Note

In contrast to the classic COD Cuvette Test (COD classic) the HT-COD is characterised by a higher digestion temperature and shorter digestion time. Users are advised to carry out a comparison with the COD classic, in order to be sure that the results obtained from their own samples when using the HT-COD are comparable to the standard.

Special note

For optimal stability until it's expiry date, it is recommended the reagent LCK 514 is stored in a fridge.

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Datatabel · Data table · Veri tablosu

DR 2800 / DR 3800 / DR 3900 / DR 5000 / DR 6000 06/2013

Software Download: www.hach-lange.com

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EN

1. Bring the sediment into suspension by inverting a few times.
2. **Carefully** pipette 2.0 mL sample.
3. Close cuvette, thoroughly clean the outside.
4. Invert.
5. Heat in the thermostat.
HT 200 S: in standard program **HT for 15 min COD classic: 2 h at 148 °C**
6. Remove the **hot** cuvette.
a. **HT 200 S:** After the lock opens, **carefully** invert **twice**.
b. **COD classic:** **Carefully** invert **twice**.
7. Allow to cool to room temperature.
a. **HT 200 S:** in the thermostat
b. **COD classic:** in a cooling rack
8. **HT 200 S:** Sediment must be completely settled before evaluation is carried out. Clean the outside of the cuvette and evaluate.
COD classic: Clean the outside of the cuvette and evaluate.

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EN: If any of the following barcode instruments is used, an automatic evaluation is carried out after the sample cuvette is inserted:

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